

360D Barbed COG



We developed the 360D COG to have **barbs more than other cogs** from all sides and **bidirectional** to give more contact and much stronger tach to the tissue, resulting in **stronger lifting**.

Thanks to our researches and development vision, The perfect designing is not just about number of barbs, but to reach;

The maximum **barbs number**.


Keeping the maximum **tensile strength**.

For maximum **lifting result and duration**.



XLINE 360D INFINITY COG

The barbs are:
Larger
Stronger
Wider



OTHER BRANDS 6D BARBED COG

The barbs are:
Smaller
Weaker
Narrow



Cross Section of XLine 360D Cog


Why Xline threads are stronger? and more durable?

Not all barbed cogs are the same;

We get a Polyfilament Barbed Cog Threads, by compressing & linking few monofilaments together.

The more number of monofilaments, and stronger linkage in between, the more strength & more duration we achieve.

Barbed Cogs extracted after 60 days



Xline 360D Cog



Other brand Cog